


# MAKOplasty UKA Results in Less Pain, More Accurate Alignment Than Those Receiving Manually Placed Oxford® Implants<sup>1</sup>

A Randomized Controlled Trial Comparing Accuracy and Outcomes of Robotic Arm Assisted UKA to Manual UKA

Investigators: Mark J Blyth, Julie Smith, Bryn Jones, Angus D MacLean III, Iain Anthony, Phillip Rowe – University of Strathclyde and Glasgow Royal Infirmary



X-Rays from first nine consecutive cases for each group

## Study Summary

- Initial results of first 100 patients from the prospective, single center, randomized controlled trial (RCT) comparing MAKOplasty UKA (RESTORIS® MCK implants) to manually placed Biomet Oxford\* implants
- VAS (Visual Analog Scale) was used to measure patient-reported pain levels
- Alignment outcomes were based on a post-operative CT protocol, which compared planned versus achieved alignment in three different planes for both the tibia and femur


## Results

- MAKOplasty resulted in significantly lower post-operative pain from day 1 to 8 weeks post-op ( $p < 0.05$ )
- MAKOplasty resulted in higher accuracy in all dimensions measured, with statistical differences ( $p < 0.01$ ) in four of the six dimensions
- A significantly higher percentage of MAKOplasty patients had excellent American Knee Society Scores at three months ( $p = 0.0031$ )


\*Oxford® is a registered trademark of Biomet, Inc.

# MAKOplasty® UKA Has Superior Alignment Compared To Manually Placed Oxford® Implants<sup>1</sup>


### Femur Varus/Valgus




### Femur Flexion/Extension



### Tibia Internal/External Rotation



### Tibial Slope



- Alignment assessment included six different parameters:
  - Femoral and tibial varus/valgus alignment, femoral flexion and tibial slope alignment, and femoral and tibial rotation
- Errors in tibial component slope have been shown to be predictors of implant failure<sup>2</sup>
- MAKOplasty UKA cases were compared to surgical plan
- Manual Oxford cases were compared to target values in surgical technique

\*Oxford® is a registered trademark of Biomet, Inc.

#### Reference

1. Blyth MJ, Smith J, Jones B, MacLean III AB, Anthony, Rowe P. Does robotic surgical assistance improve the accuracy of implant placement in unicompartmental knee arthroplasty? AAOS 2013 Annual Meeting, Chicago, IL.
2. Hernigou P, Deschamps G. Posterior slope of the tibial implant and the outcome of unicompartmental knee arthroplasty. *J Bone Joint Surg Am.* 2004;86(3):506-511.

All claims of product performance and indications for use contained within this document relate only to data submitted to and reviewed by regulatory authorities in those jurisdictions in which clearance(s) and/or approval(s) have been obtained, including the United States. No product performance claims or indications for use are made for jurisdictions in which such clearance(s) and/or approval(s) have not been obtained.



Restoring Quality of Life Through Innovation®

2555 Davie Road | Fort Lauderdale, FL 33317 | 866.647.6256 | makosurgical.com

© 2013 MAKO Surgical Corp. 208634 r00 03/13